

REMIREMONT



Remiremont is a tubular cowl with a lace pattern that creates delicate, flowing lines. Knitted with a strand of merino held together with a strand of silk mohair, it combines softness, loft and lightness.

The cowl is worked in the round from a provisional cast-on. The lace pattern repeats until the desired length is reached. After blocking, the cowl is twisted once, then the two ends are grafted together to form a continuous loop.

A cosy piece to wrap around your neck, bringing together the warmth of the yarns and the delicacy of lace.

#margueenee

#remiremontcowl



MarGueNée · Remiremont

 Instructions in English

Your pattern, step by step

Scan to discover the interactive version and pick up where you left off.

● margueenee.com

SIZE

One size.

Circumference: 35½ in ; 90 cm.

Height: 13¾ in ; 35 cm.

Adjust the cowl measurements to your liking. Keep in mind that any modifications will affect the yardage/meterage listed.

YARN

Knitting For Olive "Merino" + "Soft Silk Mohair".

3 balls of Merino + 3 balls of Mohair.

Sample shown in colorway "Iced Blue".

Or approx. 710 yd ; 650 m of each yarn at matching gauge.

GAUGE & NEEDLES

18 sts / 26 rows = 4 in ; 10 cm.

Lace pattern, washed and blocked.

Needles: US 8 ; 5 mm.

Adjust needle size to achieve correct gauge.

NOTIONS

Stitch marker.

A second US 8 ; 5 mm needle for grafting.

Crochet hook and waste yarn for the provisional cast on.

Tapestry needle.

US 8 ; 5 mm circular needle, 24 in ; 60 cm cable.

GLOSSARY

k2tog: knit the next 2 sts together.

k: knit.

p: purl.

ssk: slip slip knit.

st(s): stitch(es).

rep: repeat.

yo: yarn over.

KEYS OF THE CHART



knit



purl



yo



k2tog



ssk

TECHNIQUES

Provisional cast on:

<https://youtu.be/AKErReecbqk>

Grafting:


<https://youtu.be/dsBbmdqlaHU>

INSTRUCTIONS

Cast on **128 sts** (8 pattern repeats) with a crochet provisional cast on, **place marker** and join in the round. Repeat the **32 rounds** of the lace pattern until the cowl measures 35½ in ; 90 cm from the cast-on edge, ending with **round 32**.

Place the live stitches on holder. Cut yarn. Leave a length of yarn that measures approximately 86½ in ; 220 cm (4 times the total circumference of the provisional cast on) for the grafting.

Wet block the cowl and lay flat to dry to achieve the correct dimensions. Once completely dry move live stitches back onto a needle from waste yarn. Remove provisional cast on and arrange the stitches on an additional needle. Twist the cowl once, then graft the cast-on edge to the last round using the Kitchener stitch, following the lace pattern.



LACE PATTERN

R1: rep *p4, **k2tog**, k2, yo, p1, yo, k2, **ssk**, p3 * to end.

R2: rep *p4, k4, p1, k4, p3 * to end.

R3: rep *p3, **k2tog**, k2, yo, k1, p1, k1, yo, k2, **ssk**, p2 * to end.

R4: rep *p3, k5, p1, k5, p2 * to end.

R5: rep *p2, **k2tog**, k2, yo, k2, p1, k2, yo, k2, **ssk**, p1 * to end.

R6: rep *p2, k6, p1, k6, p1 * to end.

R7: rep *p1, **k2tog**, k2, yo, k3, p1, k3, yo, k2, **ssk** * to end.

R8: rep *p1, k7, p1, k7 * to end.

R9: rep *p1, k3, **k2tog**, k2, yo, p1, yo, k2, **ssk**, k3 * to end.

R10: rep *p1, k6, p3, k6 * to end.

R11: rep *p1, k2, **k2tog**, k2, yo, p3, yo, k2, **ssk**, k2 * to end.

R12: rep *p1, k5, p5, k5 * to end.

R13: rep *p1, k1, **k2tog**, k2, yo, p5, yo, k2, **ssk**, k1 * to end.

R14: rep *p1, k4, p7, k4 * to end.

R15: rep *p1, **k2tog**, k2, yo, p7, yo, k2, **ssk** * to end.

R16: rep *p1, k3, p9, k3 * to end.

R17: rep *p1, yo, k2, **ssk**, p7, **k2tog**, k2, yo * to end.

R18: rep *p1, k4, p7, k4 * to end.

R19: rep *p1, k1, yo, k2, **ssk**, p5, **k2tog**, k2, yo, k1 * to end.

R20: rep *p1, k5, p5, k5 * to end.

R21: rep *p1, k2, yo, k2, **ssk**, p3, **k2tog**, k2, yo, k2 * to end.

R22: rep *p1, k6, p3, k6 * to end.

R23: rep *p1, k3, yo, k2, **ssk**, p1, **k2tog**, k2, yo, k3 * to end.

R24: rep *p1, k7, p1, k7 * to end.

R25: rep *p1, yo, k2, **ssk**, k3, p1, k3, **k2tog**, k2, yo * to end.

R26: rep *p2, k6, p1, k6, p1 * to end.

R27: rep *p2, yo, k2, **ssk**, k2, p1, k2, **k2tog**, k2, yo, p1 * to end.

R28: rep *p3, k5, p1, k5, p2 * to end.

R29: rep *p3, yo, k2, **ssk**, k1, p1, k1, **k2tog**, k2, yo, p2 * to end.

R30: rep *p4, k4, p1, k4, p3 * to end.

R31: rep *p4, yo, k2, **ssk**, p1, **k2tog**, k2, yo, p3 * to end.

R32: rep *p5, k3, p1, k3, p4 * to end.